

SURFACE / READY TO USE DISINFECTANT

## DESCOCEPT SENSITIVE



PRESERVE

ALCOHOLIC RAPID DISINFECTANT

1 / 10

Efficient disinfection and cleaning performance  
**with long service life**

Wide range of applications  
**for daily work**

Very good material compatibility  
**with low ethanol content**

Ideally suited  
**for use in the food industry**



Use in HACCP  
concepts

SURFACE / READY TO USE DISINFECTANT

DESCOCEPT  
SENSITIVE



PRESERVE

PRODUCT DESCRIPTION

2 / 10

## Versatile use on surfaces and medical devices

**DESCOCEPT SENSITIVE** is a ready-to-use, alcohol-based rapid disinfectant based on ethanol for gentle disinfection and cleaning of medical devices and medical equipment as well as surfaces of all types. The product is particularly suitable for areas near patients in which a broad spectrum of activity is required.

Besides providing effective disinfection within a very short time, the product has a very good material compatibility.

### APPLICATIONS AND NOTES

#### According to Biocidal Products Regulation (BPR)

Rapid disinfection and cleaning of alcohol-resistant medical inventory and all types of surfaces, mainly in near patient areas.

#### According to EU Medical Devices Regulation

For rapid disinfection and cleaning of non invasive medical devices, especially suited for sensitive surfaces of medical devices.

#### Application

To ensure complete wetting apply the undiluted solution evenly on alcohol-resistant surfaces and keep them moist for the required contact time.

If the surface is highly contaminated, the material should be removed mechanically immediately and surface disinfection should only be carried out afterwards.

According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

#### Application notes

DESCOCEPT SENSITIVE can be used for wipe disinfection or as a spray. Suitable for use in the food sector.

#### Composition

100 g solution contain: 45 g ethanol.

#### Material compatibility

Wide range of applications on surfaces and medical devices. (see page 5 - 8)

#### Precautionary and hazard statements

Flammable liquid and vapour. If medical advice is needed, have product container or label at hand. Keep out of reach of children. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Do not breathe vapour/spray. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents/container to approved disposal company or local collection.

For professional use only by personnel with corresponding specialist knowledge according to national directives.

**Use disinfectants safely. Always read label and product information before use.**

## SURFACE / READY TO USE DISINFECTANT

# DESCOSEPT SENSITIVE



PRESERVE

### SPECTRUM OF EFFICACY AND CONTACT TIMES

3 / 10



SPECTRUM OF EFFICACY AND CONTACT TIMES			30 s	1 min	2 min	3 min
<b>application recommendations for surface disinfection (wiping disinfection)</b>						
bactericidal, yeasticida <sup>1)</sup>	VAH/EN <sup>2,3</sup>	with mechanical action, dirty conditions			•	
tuberculocidal (M. terrae)	EN 14348	clean and dirty conditions		•		
mycobactericidal (M.terrae, M.avium)	EN 14348	clean and dirty conditions			•	
limited spectrum virucidal	EN 14476	clean conditions				•
active against enveloped viruses	EN 14476	clean and dirty conditions	•			
<b>application recommendations for surface disinfection (spraying disinfection)</b>						
bactericidal, yeasticidal <sup>1,4</sup>	EN 17387	with mechanical action, dirty conditions		•		
bactericidal, yeasticidal <sup>1,5</sup>	EN 13697	without mechanical action, clean and dirty conditions		•		
tuberculocidal (M. terrae)	EN 14348	clean and dirty conditions		•		
mycobactericidal (M.terrae, M.avium)	EN 14348	clean and dirty conditions			•	
limited spectrum virucidal	EN 14476	without mechanical action, clean conditions				•
	EN 16777	without mechanical action, clean conditions				•
active against enveloped viruses	EN 14476	without mechanical action, clean conditions		•		
	EN 16777	without mechanical action, clean conditions		•		

1 – including phase 2 step 1 and phase 2 step 2 tests (quantitative suspension tests and quantitative carrier tests)

2 – EN 16615, EN 13727, EN 13624, VAH Methode 8

3 – VAH listed 5 minutes

4 – medical

5 – industry, food, household, public institutions

SURFACE / READY TO USE DISINFECTANT

## DESCOCEPT SENSITIVE



PRESERVE

### SPECTRUM OF EFFICACY AND CONTACT TIMES

4 / 10



SPECTRUM OF EFFICACY AND CONTACT TIMES			30 s	1 min	2 min	3 min
<b>additional results</b>						
bactericidal	EN 16615	dirty conditions			•	
	EN 13697	clean conditions	•			
		dirty conditions		•		
	EN 13727	dirty conditions		•		
yeastical (C. albicans)	EN 1276	dirty conditions	•			
	EN 16615	dirty conditions			•	
		clean and dirty conditions	•			
	EN 13697	clean and dirty conditions		•		
		dirty conditions		•		
active against adenoviruses	EN 1650	dirty conditions	•			
	EN 14476	clean conditions				•
	EN 16777	without mechanic action, clean conditions				•
active against noroviruses (MNV)	EN 14476	dirty and clean conditions		•		
	EN 16777	without mechanic action, dirty conditions				•
active against vaccinia viruses	EN 14476	dirty and clean conditions	•			
	EN 16777	without mechanic action, clean conditions		•		
active against corona viruses	EN 14476	dirty conditions	•			

SURFACE / READY TO USE DISINFECTANT

# DESCOCEPT SENSITIVE



PRESERVE

## MATERIAL COMPATIBILITY

5 / 10



MATERIAL METALS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
Stainless steel V2A			<ul style="list-style-type: none"> <li>Medical transport chairs</li> <li>Toilet chairs</li> <li>Rollators</li> <li>Walking frames</li> </ul>	
Copper			<ul style="list-style-type: none"> <li></li> </ul>	
Anodized aluminum			<ul style="list-style-type: none"> <li>Instrument containers</li> <li>Sterile containers</li> </ul>	
Brass		<ul style="list-style-type: none"> <li></li> </ul>	<ul style="list-style-type: none"> <li>Luer lock adapter</li> <li>Stethoscope (chrome-plated brass)</li> </ul>	
Aluminum	<ul style="list-style-type: none"> <li></li> </ul>			

SURFACE / READY TO USE DISINFECTANT

# DESCOCEPT SENSITIVE



PRESERVE

## MATERIAL COMPATIBILITY

6 / 10



MATERIAL PLASTICS: ELASTOMERS	not recommended	limited recommended	recommended	APPLICATION PRODUCT EXAMPLE
CR (chloroprene rubber)			•	
PUR (polyurethane)			•	Medical transport chairs
EPDM (ethylene propylene diene (monomer) rubber)			•	Nursing trolleys
MQ (polydimethylsilicone rubber)			•	
VMQ (polymethylsiloxane vinyl rubber)			•	Respiratory masks
PVMQ (polymethylsiloxane-phenyl vinyl rubber)			•	Medical keyboards and mice

SURFACE / READY TO USE DISINFECTANT

# DESCOCEPT SENSITIVE



PRESERVE

## MATERIAL COMPATIBILITY

7 / 10



MATERIAL PLASTICS: THERMOPLASTICS	not recommended	limited recommend	recommended	APPLICATION PRODUCT EXAMPLE
PC (polycarbonate e.g. Makrolon)			<ul style="list-style-type: none"> <li>• Ultrasound devices</li> <li>• ECG devices</li> <li>• EEG devices</li> <li>• MRI devices</li> <li>• CT devices</li> </ul>	
PC/ABS (polycarbonate/acrylonitril-butadiene-styrene)			<ul style="list-style-type: none"> <li>• X-ray devices</li> <li>• HF surgery housing components</li> <li>• Abdominal ultrasound probes</li> <li>• Perfusors</li> <li>• Incubators</li> </ul>	
ABS (acrylonitril-butadiene-styrene)			<ul style="list-style-type: none"> <li>• Patient monitoring monitors</li> <li>• Medical keyboards and mice</li> </ul>	
PE-HD (polyethylene-High Density)			<ul style="list-style-type: none"> <li>• Storage and transport containers</li> </ul>	
POM (polyoxymethylene)			<ul style="list-style-type: none"> <li>• Seat shells</li> </ul>	
PA (polyamide)			<ul style="list-style-type: none"> <li>• X-ray chute</li> </ul>	
EPP (expanded (foamed) polypropylene)			<ul style="list-style-type: none"> <li>• Fascia rollers</li> </ul>	
PET (polyethylene terephthalate)		•		
PMMA (polymethylmethacrylate)	•			Acrylic and plexiglass, incubators



SURFACE / READY TO USE DISINFECTANT

# DESCOCEPT SENSITIVE



PRESERVE

## MATERIAL COMPATIBILITY

8 / 10



MATERIAL TEXTILES	not recommended	limited recommend	recommended	APPLICATION PRODUCT EXAMPLE
synthetic leather (e.g. with PVC and PUR top layer)				<ul style="list-style-type: none"> <li>Examination and treatment couches</li> <li>Storage aid</li> </ul>
				<ul style="list-style-type: none"> <li>Gynaecological examination chair</li> <li>Dental chairs</li> </ul>
				<ul style="list-style-type: none"> <li>Operating tables</li> <li>Patient transport couches</li> </ul>



SURFACE / READY TO USE DISINFECTANT

## DESCOCEPT SENSITIVE



PRESERVE

### PURCHASING INFORMATION

9 / 10

Product	Single unit	Delivery	Content	Art. Nr.
	Spray Bottle	6	750 ml	00-323DS-0075-02
	Bottle	10	1 L	00-323DS-010
	Canister	3	5 L	00-323DS-050

National information may differ. For further information, please contact our subsidiary or your local dealer.  
The availability of the products and container sizes depends on a completed national registration.



### CERTIFICATIONS

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001, has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.

SURFACE / READY TO USE DISINFECTANT

## DESCOCEPT SENSITIVE



PRESERVE

### PRODUCT FAMILY OVERVIEW

10 / 10



DESCOCEPT SENSITIVE WIPES



DESCOCEPT SENSITIVE WIPES



DESCOCEPT SENSITIVE WIPES XL



DESCOCEPT SENSITIVE



DESCOCEPT SENSITIVE



### RECOMMENDED NON-WOVEN WIPE DISPENSER SYSTEMS



ONE SYSTEM+ PLUS



ONE SYSTEM BASIC



DESCO WIPES



ECO WIPES