

INFRARED FOREHEAD THERMOMETER

Non-contact Measurement "O" Risk of Cross-infection



LFR30B




SAFE



FAST



CONVENIENT



SPECIFICATIONS

Measuring Position	Forehead
Measuring Distance	Ocm~5cm
Display	LCD Screen
Resolution	0.1°C / 0.1°F
Measuring Range	Human Body: 32.0°C ~ 43.0°C Object: 0°C ~ 100°C
Accuracy	For body: ±0.3°C (34°C~43°C) ±0.4°C (32°C~33.9°C)
Operating Temperature & Humidity	Temperature: 5°C~40°C Humidity: ≤85%RH
Storage Temperature & Humidity	Temperature: -20°C~55°C Humidity: ≤95%RH
Power Supply	2 X 1.5V AAA bateries
Service Life	5 years

Lepu Medical Technology (Beijing) Co., Ltd. -