


HYGIENE PRODUCTS MANUFACTURE





Bournas Medicals has been operating for more than **30 years** in the field of medical disposable products. In 2006, the company took a step forward starting the **production** of hygiene paper products for the medical, dental, beauty and cleaning sector.

Bournas Medicals holds **high market share in Greece** by covering the entire range of health products with distributors all around the country. From 2010, the company is **exporting** its products worldwide by having significant presence in more than **40 countries** like France, Belgium, Germany, Australia, Cyprus, Lebanon, Syria, Qatar etc.

Since the foundation in 1988, the company's vision has always been the focus on **investing** its resources on the development of new forms of **manufacturing** technics thus providing products of **excellence**. This is accomplished by focusing, always, on customer's needs and requirements.

Trust is the philosophy of our business. Bournas Medicals' **goal** is to provide excellent **pre and after-sales** service by its highly **experienced personnel**, offering assistance wherever it is required.

The company is **certified** with ISO 9001:2015 for production, trading and distribution of all products. With special emphasis and attention to detail, we look after the quality of our products so that we can meet the strictest specifications and controls required by European Certification (CE).

Company's main **product categories** include hygiene examination rolls & tissue paper products, dental rolls & bibs, airlaid, sanitary & incontinence products, disposable gloves, cardiology, protection, hospital & laboratory products, among other medical disposables.

The **wide range** of products cover many **sectors** in terms of **usage**, such as hospitals, clinics, private doctors, dentists, beauty salons, spa & massage centers, orthopedics, maternity & gynecological clinics, laboratories & diagnostic centers, tattoo studios, cleaning companies, food industry, etc.

Bournas Medicals is highly committed to its customers on offering finest quality products. First and foremost important concern, is to provide full satisfaction by offering the best possible quality products in the most competitive prices.



ELITE
Thermal Laminated Rolls (100% pure pulp)

PREMIUM
Laminated Rolls with Glue (100% pure pulp)

ECO
Laminated Rolls with Glue (100% recycled pulp)

CLASSIC
Plain Paper Rolls (100% pure pulp)

REGULAR
Paper Rolls, embossed (100% pure pulp)

SUPREME
Non Woven Rolls (100% spunbond)

PRESTIGE
Non Woven + PE Laminated Rolls (100% spunbond)

GENTLE
Tissue paper wipes/ Folded or un folded

HYPER
Airlaid Rolls/embossed or Not & Airlaid wipes/ Folded

DELICATE
Towel Industrial Rolls/ Embossed or Not

RAINBOW
Dental Rolls, Dental Head Rest Covers,
Dental Tray Paper & Dental Bibs (100% pure pulp)



ELITE

Thermal Laminated Rolls (100% pure pulp + PE)

Thermal laminated examination roll of 100% pure pulp, waterproof, absorbent, high resistant, perforated for easy cutting. This roll is manufactured out of a PE film placed at the bottom of the roll and a sheet of paper of exceptional quality and strength, placed at the top of the roll. The 2 plies stick together with thermal lamination process which makes the roll stronger and the 2 plies very hard to separate.

All rolls are individually shrinked for safety reasons.



Usage

hospitals, clinics,
private doctors,
beauty salons,
generally wherever
liquids are used



white

green

light blue

pink

purple

yellow

black

print

50cm	50m							
	80m							
	100m							
58cm	50m							
	80m							
	100m							
68cm	50m							
	80m							
	100m							

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.





PREMIUM

Laminated Rolls with Glue

(100% pure pulp + PE)

Glue laminated examination roll of 100% pure pulp, waterproof, absorbent, high resistant, perforated for easy cutting. This roll is manufactured of a PE film placed at the bottom of the roll and a sheet of paper of exceptional quality and durability placed on the top side of the roll. This examination roll prevents possible leaks during the use of liquids and also offers a sense of comfort and cleanliness.

All rolls are individually shrinked for safety reasons.

Usage

hospitals, clinics,
private doctors,
beauty salons,
generally wherever
liquids are used



white

green

light blue

pink

purple

yellow

black

print

50cm	50m	PRINT	2	38cm	132	12 or 6	32 or 64
	80m	PRINT	2	38cm	210	8	32
	100m	PRINT	2	38cm	264	6	32
58cm	50m	PRINT	2	38cm	132	12 or 6	32 or 64
	80m	PRINT	2	38cm	210	8	32
	100m	PRINT	2	38cm	264	6	32
68cm	50m	PRINT	2	38cm	132	12	24
	80m	PRINT	2	38cm	210	8	24
	100m	PRINT	2	38cm	264	6	24

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.





ECO

Laminated Rolls with Glue (100% recycled pulp + PE)

Glue laminated examination roll from recycled pulp of excellent quality and PE coating for waterproof protection. With a perforated surface for easy cutting, this roll consists of a PE film placed on the bottom of the roll and a sheet of recycled paper on the top side, offering durability and high protection.

All rolls are individually shrinked for safety reasons.



Usage

public hospitals,
clinics,
generally wherever
liquids are used


white

Roll Width	Roll Length	Perforations	Rolls per Box	Box Dimensions (cm)	Rolls per Pallet	Pallet Dimensions (cm)	
50cm	50m	<input type="checkbox"/> PRINT	2	38cm	132	12 or 6	32 or 64
	80m	<input type="checkbox"/> PRINT	2	38cm	210	8	32
	100m	<input type="checkbox"/> PRINT	2	38cm	264	6	32
60cm	50m	<input type="checkbox"/> PRINT	2	38cm	132	12 or 6	32 or 64
	80m	<input type="checkbox"/> PRINT	2	38cm	210	8	32
	100m	<input type="checkbox"/> PRINT	2	38cm	264	6	32
70cm	50m	<input type="checkbox"/> PRINT	2	38cm	132	12	24
	80m	<input type="checkbox"/> PRINT	2	38cm	210	8	24
	100m	<input type="checkbox"/> PRINT	2	38cm	264	6	24

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.





CLASSIC & REGULAR

Paper Rolls plain or embossed (100% pure pulp)

Thick and absorbent 2ply paper examination roll, manufactured by first class quality pulp to protect medical examination couch as well as for any other professional and domestic use. This roll is soft, absorbent and it can be easily cut on a customized perforation. This roll is available in Plain (Classic) and Embossed (Regular) paper.

All rolls are individually shrinked for safety reasons.



Usage

hospitals,
clinics,
private doctors,
beauty salons



white

light blue

		PLAIN		EMBOSSED							
Width	Length	PRINT	PRINT	PRINT	PRINT	Perforation	Rolls per box	Rolls per pallet	Rolls per pallet		
50cm	50m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	132	12 or 6	32 or 64
	80m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	210	8	32
	100m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	264	6	32
58cm	50m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	132	12 or 6	32 or 64
	80m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	210	8	32
	100m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	264	6	32
68cm	50m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	132	12	24
	80m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	210	8	24
	100m	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	PRINT	2	38cm	264	6	24
50cm	50m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	132	12 or 6	32 or 64
	80m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	210	8	32
	100m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	264	6	32
58cm	50m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	132	12 or 6	32 or 64
	80m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	210	8	32
	100m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	264	6	32
68cm	50m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	132	12	24
	80m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	210	8	24
	100m	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	PRINT	1	38cm	264	6	24

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.





SUPREME & PRESTIGE

Non-Woven (TNT) Rolls & Non-Woven Laminated Rolls (100% spunbond) & (100% spunbond + PE)


Supreme non-woven Examination Roll is made of 100% polypropylene, hydrophobic, non-woven material. It is perforated for easy cutting. Suitable for testing beds or beds therapy in clinics, physiotherapy, massage and spa centers.

Prestige non-woven & PE 2ply Examination Roll is thick and absorbent. This Roll is soft yet of high resistance, suitable to cover the medical examination couch. It is perforated for easy cutting. Ideal for medical, clinics, full body therapy centers.



Usage

hospitals, clinics,
private doctors, beauty
salons, spa, massage
& therapy centers



	Weight	Width	Length	Perforation	Rolls per box	Rolls per pallet
15gr	50cm	70m	1	38cm	184	12
	58cm	70m	1	38cm	184	12
	68cm	70m	1	38cm	184	12
20gr	50cm	70m	1	38cm	184	9
	58cm	70m	1	38cm	184	9
	68cm	70m	1	38cm	184	9
40gr	50cm	70m	1	38cm	184	6
	58cm	70m	1	38cm	184	6
	68cm	70m	1	38cm	184	6
NON-WOVEN & PE	50cm	50m	2	38cm	132	15
	58cm	50m	2	38cm	132	15
	68cm	50m	2	38cm	132	15

Custom width & length upon request.

Custom perforation upon request.





HYPHER

Airlaid Rolls & Wipes

Air Laid cleaning Rolls & Wipes can absorb 7 times their weight. They can be rinsed and reused up to 9 times. They are soft, absorbent, of high resistance and can be sterilized. Airlaid products are non-linting and strong, even when wet. They have a textile-like surface and drape.







Usage

hospitals, private doctors, beauty salons, food sector, industrial sector, hygiene & cleaning industry, domestic use, personal care

	Weight	Width	Length	Perforation	Rolls per box	Rolls per pallet	Rolls per pallet	Rolls per pallet	
AIRLAID WIPES	50gr	30cm	33cm	PRINT	1	50	24	32	
	50gr	30cm	40cm	PRINT	1	50	24	32	
	50gr	40cm	40cm	PRINT	1	50	18	32	
	50gr	40cm	70m	PRINT	1	50	12	32	
	80gr	40cm	70m	PRINT	1	50	6	32	
AIRLAID EXAMINATION ROLLS	50cm	25m		PRINT	1	38cm	132	12 or 6	32 or 64
	58cm	25m		PRINT	1	38cm	132	12 or 6	32 or 64
AIRLAID INDUSTRIAL ROLLS	500gr	25cm	40m	PRINT	1	38cm	105	12	-
	1000gr	29cm	68m	PRINT	1	38cm	179	6	-
	2000gr	29cm	138m	PRINT	1	38cm	363	2	-

Custom width & length upon request.

Custom perforation upon request.






GENTLE

Tissue Paper Sheets

Tissue Paper sheets are made from 100% pure or recycled pulp. They are soft and durable in use, liquid absorbent and are packed in successive dimensional sheets of 40cm x 60cm. These paper sheets leave no residue and do not irritate the skin. Available in packs of 5 kg, 1 kg and 300 grams in durable plastic bag. Folded inside strong carton box for complete safety during transport and storage, also available.



Usage

hospitals, clinics,
private doctors,
beauty salons, all
professional sectors

ALFA - ALFA	100% pure pulp, 86% brightness (±5% tolerance)	40cm	60cm	5kgr	4	80 PACKAGES
		40cm	60cm	1kgr	10	32
		40cm	60cm	300gr	25	40
		40cm	60cm	5 kgr folded	1	56
EXTRA - SOFT	recycled pulp, 67% brightness (±5% tolerance)	40cm	60cm	5kgr	4	80 PACKAGES
		40cm	60cm	1kgr	10	32
		40cm	60cm	300gr	25	40
		40cm	60cm	5 kgr folded	1	56
ALFA - SOFT	recycled pulp, 55% brightness (±5% tolerance)	40cm	60cm	5kgr	4	80 PACKAGES
		40cm	60cm	1kgr	10	32
		40cm	60cm	300gr	25	40
		40cm	60cm	5 kgr folded	1	56





DELICATE

Towel Industrial Rolls & Wipes

Towel Industrial Rolls (Plain or Embossed). The use of towel is a fast, effective, convenient and economical way of improving hand hygiene. The extended range of hand wiping products includes larger, more absorbent towels than usual, making the hand drying not only a pleasant experience, but also quick and effective.



Usage
hospitals, clinics, private doctors, beauty salons, all professional sectors


zig zag wipes	21cm	25cm	white	1	4000	32

500gr roll	25cm	52m	white	2	38cm	137	12
1000gr roll	29cm	104m	white	2	38cm	274	6
4000gr roll	29cm	400m	white	2	38cm	1053	2

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.





RAINBOW

Dental Rolls with or w/o neckline (100% pure pulp + PE)

Dental Bibs on Rolls with perforation for easy cutting. They are waterproof, soft, of high resistance, strong, made of 100% pure pulp and a sheet of PE. Dental wipes are also available with neckline (perforated).

All rolls are individually shrinked for safety reasons.



Usage

dentists, dental clinics, hospitals, private doctors, dental usage mostly

DENTAL BIBS ON ROLL WITH NECKLINE



53cm

60cm




2

80pcs per roll

9

48

DENTAL BIBS ON ROLL



29cm

50m



2

38cm

132

24 or 12

32 or 64

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.





RAINBOW

Dental Head Rest Covers

Dental Tray Paper

Head Rest Covers are made of 1ply paper and 1ply PE. Used for covering the head rest area of the dentist's chair. Dental usage mostly.


Dental Tray paper is made of high absorbent paper. Used for placing the dentist's tools on the tray. Dental usage mostly.



Usage

medical & dental clinics, private doctors, dental usage mostly

HEAD REST COVERS



29cm

30cm




2

150pcs

8 x 150 pcs

20

DENTAL TRAY PAPER



18cm

28cm



1

250pcs

10 x 250 pcs

80





RAINBOW

Dental Bibs (100% pure pulp + PE)

Dental Bibs made of 100% pure paper pulp & PE.
Available in 2 plies & 3 plies:

- 2plies: 1 sheet of paper & 1 sheet of PE
- 3plies: 2 sheets of paper & 1 sheet of PE

Dental Bibs are soft, waterproof, absorbent, of high resistance.



Usage

dentists, hospitals, clinics, private doctors, beauty salons, manicure, pedicure, physiotherapy, massage & spa centers



white green light blue pink yellow purple lime bordeaux blue orange fuchsia green black print

33cm	45cm		3 2 ply Paper + 1 ply PE	inner dispenser	125pcs	500pcs	72	
				PE bag	500pcs	500pcs	72	
				2 1 ply Paper + 1 ply PE	inner dispenser	125pcs	500pcs	72
					PE bag	500pcs	500pcs	72

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.



Dental Bibs are packed in cartons of 500pcs.

There are two options of packaging in cartons:

- Either protected by a PE bag
- or in paper dispensers per 125pcs.

Dispensers are paper boxes perforated in such a sophisticated way that can be opened half way and/or all the way down in a delicate curved shape. The opening suits for easy disposal, while the half way opening option is ideal for safety and protection reasons.





BOURNAS MEDICALS:

📍 Seneka 18, Kifisia 145 64 Athens, Greece

☎ +30 210 2622 568

☎ +30 210 2612 089

✉ info@bournas-medicals.gr

🏠 www.bournas-medicals.gr

📷 [@bournas.medicals](https://www.instagram.com/bournas.medicals)

📘 [@bournas.medicals](https://www.facebook.com/bournas.medicals)