



**BOURNAS**  
MEDICALS

Health supplies since 1988

Production • Imports • Exports

# PRODUCT CATALOGUE

SOFTcare

 **dental** series




**SOFTcare**  
dental



A complete series  
of dental disposable products  
guarantees healthy smiles

Certified Protection



[www.bournas-medicals.gr](http://www.bournas-medicals.gr)

Dental Bibs Standard & Plus 4

Dental Bibs on a Roll 6

Dental Bibs on a Roll with Neckline 7

Dental Head Rest Covers 8

Dental Tray Papers 9

Dental Saliva Rejectors 10

Dental Cups 10

Dental Cotton Rolls 11

Comparison Chart: Bibs & Rolls 14



## ABOUT US:

“Bournas Medicals” is a modern **Greek** company, operating since 1988 in the production and trade of **medical disposable products** and **health supplies**. A **leader** company with **productive** and **import** activity, while it has a dynamic presence with active **export activity** in more than **50 countries worldwide**.

With our extensive catalogue of certified products, we provide complete solutions covering the needs of every business or corporate entity that trades and distributes medical and health consumables. Thus, indirectly but effectively we support the services of various sectors:

- Medical and hospital care
- Health and safety
- Beauty and wellness.

### Our Vision:

We **inspire** companies related to our products all **around the world**. We are the **No. 1 choice** for the wholesaler of health supplies.

### Our Mission:

We **support** the modern wholesale **business** as **one stop point** for the **supplies of medical and health disposable products**.

### Our Goals:

- To provide **products** that **meet** the **strictest criteria** and **certifications**.
- To **adjust our pricing policy** and **catalogue** in a way **that ensures the increase of sales for our customers**.
- The creation of **certified affordable** medical consumables, through our own brand label (“**SOFT Care**”), based on the needs and trends of the markets we operate.
- The continuous **expansion** of our **catalogue** with new products & options.
- To support our **B2B** network of **customers worldwide**, with high level after-sales services and training for our products.

## OUR PRODUCTION LINE AT YOUR SERVICE

Standard and bespoke products for the modern wholesale business.





## DENTAL

### Bibs Standard & Plus

- **2ply or 3ply:** 100% pure paper pulp & HDPE
- Glue laminated
- Super soft & highly resilient
- Water proof
- Extra wet strength
- Folded
- Packaging convenient dispensers or in PE bag

CE MD

#### Suggested use

DENTAL

BEAUTY

MEDICAL

TATTOO



white green light blue yellow purple lime green blue orange pink bordeaux fuchsia dark green black kids

	HEIGHT	WIDTH	COLOURS	PLIES	INNER PACKAGING	PIECES	PIECES/CARTON	CARTONS/PALLET
PLUS	33 cm	45 cm	 PRINT kids	2ply paper & PE	Dispenser	4 x 125 pcs	500 pcs	72
					PE bag	500 pcs	500 pcs	72
STANDARD	33 cm	45 cm	 PRINT kids	1ply paper & PE	Dispenser	4 x 125 pcs	500 pcs	72
					PE bag	500 pcs	500 pcs	72

Ability to print other colors or designs (e.g. logo) upon request.

## Bibs Standard & Plus Packaging

There are two options of packaging in cartons:

- A 500 pcs** protected by a PE bag inside a carton.
- B 125 pcs x 4 dispensers** inside a carton.

Our Dispensers are made of two side coated duplex board with grey back, perforated in a way that can be opened either half way and/or all the way down in a delicate curved shape.

The opening is ideal for easy dispensing, while the half way opening option helps keep the bibs better protected, yet easily accessible.



A



B



## DENTAL

### Bibs on a Roll

- **2ply:** 100% pure paper pulp & HDPE
- Glue laminated
- Super soft & highly resilient
- Water proof
- Extra wet strength



#### Suggested use



DENTAL




BEAUTY



MEDICAL



TATTOO



HEIGHT	WIDTH	COLOURS	PLIES	PIECES	ROLLS/CARTON	CARTONS/PALLET
29 cm	38 cm	PRINT	2	132pcs per roll	9	72

Custom width & length upon request.

Ability to print other colors or designs (e.g. logo) upon request.

Custom perforation upon request.






**DENTAL**

**Bibs on a Roll with Neckline**

- **2ply:** 100% pure paper pulp & HDPE
- Glue laminated
- Super soft & highly resilient
- Water proof
- Extra wet strength
- Folded and perforated



**Suggested use**



HEIGHT	WIDTH	COLOURS	PLIES	PIECES	ROLLS/CARTON	CARTONS/PALLET
53 cm	60 cm		2	80pcs per roll	9	48



## DENTAL

### Head Rest Covers

- **2ply:** 100% pure paper pulp & PE
- Thermal laminated
- Soft & highly resilient
- Water proof
- Extra wet strength
- Packaging convenient dispensers




#### Suggested use

DENTAL

BEAUTY

+ MEDICAL



HEIGHT

29 cm



WIDTH

30 cm



COLOURS



PLIES

2



PIECES

150 pcs



PIECES/CARTON

8 x 150 pcs



CARTONS/PALLET

20



Headrest Cover dispenser





**DENTAL**

**Tray Papers**

- Crepe paper
- Chlorine free
- Suitable for sterilization up to 134 °C (apart from Bordeaux)



**Suggested use**

DENTAL

BEAUTY

MEDICAL



HEIGHT	WIDTH	COLOURS	PLIES	PIECES	PIECES/CARTON	CARTONS/PALLET
18 cm	28 cm		1	250 pcs	10 x 250 pcs	80

Dental Tray paper



**DENTAL**


**Saliva Ejectors & Cups**

**Dental Saliva Ejectors**

- Made from non-toxic PVC
- Free of phthalates & latex
- Filter from LDPE
- Removable cap
- Wide variety of colours
- Maximum flexibility & stability



**Suggested use**




white	green	light blue	bordeaux	purple	yellow	fuschia	blue	orange	clear
LENGTH CANULA / FILTER	DIAMETER CANULA / FILTER	COLOURS	PIECES	PIECES/CARTON	CARTONS/PALLET				
150 / 14 mm	6,5 / 9 mm		100	50 x 100	24				

**Dental Cups**

- Made of polypropylene
- 180 ml
- Wide variety of colours



**Suggested use**



white	green	light blue	bordeaux	purple	pink	yellow	fuschia	blue	orange	black
ML	HEIGHT	DIAMETER	COLOURS	PIECES	PIECES/CARTON	CARTONS/PALLET				
180 ml	81 mm	45 mm		100	30 x 100	24				

**DENTAL**

**Cotton Rolls**

**Dental Cotton Rolls**

- Highly absorbent core
- 100% pure cotton
- Free of chlorine
- Maintain their shape
- Non-adherent
- Autoclavable



**Suggested use**



No 1 - 8mm



No 2 - 10mm



No 3 - 12mm



SIZE	LENGTH	DIAMETER	COLOURS	PIECES per 300gr (aprox)	PIECES/CARTON	CARTONS/PALLET
No 1	40 mm	8 mm	●	~900	30 x ~900	20
No 2	40 mm	10 mm	●	~680	30 x ~680	20
No 3	40 mm	12 mm	●	~450	30 x ~450	20



Partners in Health, Business & Trust

By choosing “Bournas medicals” for your medical & health supplies you have **multiple benefits**



**Experience**

More than **30 years** close to your needs









**Experience**

- Presence with our products in more than **50 countries** worldwide
- Production, import and export activity



**A variety of products**

- More than **1.400 product** codes in **15 product categories**
- Our catalogue includes the finest international brands in health supplies
- Continuous research and development
- Our products meet the needs of various sectors:

- |  |  |
|--|--|
|  Medical  |  Health |
|  Dental   |  Beauty |
|  HO.RE.CA |  Tattoo |
|  Industry |  |





## QUALITY ASSURANCE

### Certifications and standards

A number of certifications standards guarantee the quality and safety for the end user.



**Conformity marking**  
for regulating the goods sold within European Union



**Medical Devices Mark**  
All products marketed with this symbol comply with requirements stipulated in European Directives



**ISO 13485: 2016 by G-Certi**  
Certification for the Quality Management of Medical Devices



**Health Technology Certification from HTCert**  
Certification for the good practice of distribution of medical devices



**ISO 9001: 2015 by Qcert**  
Certification for the Quality Management



**ISO 14001: 2015 by Qcert**  
Certification for the Effective Environmental Management



**ISO 45001: 2018 by EBETAM MIRTEC**  
Certification for the Health and Safety at Work



**ISO 50001: 2018 by EBETAM MIRTEC**  
Certification for the Efficient Energy Management



**Registered for the use of GS1 Barcodes**  
according to the new European Regulation for Medical Devices

### Partners in Trust

- We do not sell directly to our customer's clients or end users
- We are fully in line and up to date with the existing legislations and developments in the market we operate



### No1 value for money choice

Thanks to:


- Our competitive prices
- The number of our product categories
- The vast number of our product categories
- The wide range of our certified own brand label products



















### Business to Business Client Service





We guarantee the best service thanks to:

- Our **one stop shop** character and the wide range
- Our **B2B e-shop** for even faster and easier placement of your orders
- We provide **consulting** on health supplies for the **best service of your customers**



# Comparison Chart: Bibs & Rolls


SERIES	COMPOSITION	ICONS	PLIES	LAMINATION	SOFTNESS
 <b>DENTAL BIBS PLUS</b>	100% pure paper pulp	 Comply to MDR 2017/745 Class I  	3	Glue	●●●○○
 <b>DENTAL BIBS STANDARD</b>	100% pure paper pulp	 Comply to MDR 2017/745 Class I  	2	Glue	●●●○○
 <b>DENTAL BIBS ON A ROLL</b>	100% pure paper pulp	 Comply to MDR 2017/745 Class I  	2	Glue	●●●○○
 <b>DENTAL BIBS ON A ROLL WITH NECKLINE</b>	100% pure paper pulp	 Comply to MDR 2017/745 Class I  	2	Glue	●●●○○

ABSORBENCY	RESILIENCE	PE	WATERPROOF / HYDROPHOBIC	COLOURS	SUGGESTED USE
●●●●○	●●●●○	HDPE	Waterproof		DENTAL, MEDICAL, BEAUTY, TATTOO
●●●○○	●●●○○	HDPE	Waterproof		DENTAL, MEDICAL, BEAUTY, TATTOO
●●●○○	●●●●○	HDPE	Waterproof		DENTAL, MEDICAL, BEAUTY, TATTOO
●●●○○	●●●●○	HDPE	Waterproof		DENTAL, MEDICAL, BEAUTY



Health supplies since 1988

More than **1.400 product codes** in **15 product categories**



Partners in Health, Business and Trust

18, Seneka str., Kifissia, 145 64, Athens - Greece  
tel.: +30210 2622568, +30211 2138 744  
fax: +30210 2612089

[www.bournas-medicals.gr](http://www.bournas-medicals.gr)