

• **ELITE**

Paper & PE examination rolls

- **2ply:** 100% pure paper pulp & LDPE
- Thermal laminated
- Super soft & highly resilient
- Waterproof
- Extra wet strength
- Perforated
- Wrapped separately with PE shrink film



Suggested use

- MEDICAL
- BEAUTY
- HEALTH
- TATTOO



	REF	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/CARTON	CARTONS/PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
STANDARD	108.021.ST	39 cm	50 m		2	38 cm	132	12	40	40 x 30 x 41 cm	80 x 120 x 220 cm
	108.023.ST		50 m		2	38 cm	132	12			
	108.023.STM75	50 cm	75 m		2	38 cm	210	8	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.023.STM100		100 m		2	38 cm	264	6			
	108.025.ST		50 m		2	38 cm	132	12			
	108.025.STM75	58 cm	75 m		2	38 cm	210	8	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.027.ST		50 m		2	38 cm	132	12			
	108.027.STM75	68 cm	75 m		2	38 cm	210	8	24	40 x 30 x 69 cm	80 x 120 x 225 cm



Custom width & length and perforation upon request.
For further information please contact us.

Ability to print other colors or designs
(e.g. logo) upon request.

All packaging dimensions may vary +/- 1% | All roll dimensions may vary +/- 2%
Palletization heights are suggested for transportation cost efficiency for your orders.

Certified protection

• **HYPER**

Airlaid examination rolls & wipes

“Softer than air,
strong like a textile”

- 100% airlaid
- Super soft & resilient
- Waterproof (2ply version)
- Max wet strength & absorbency
- Textile look & feel
- Non-linting
- Perforated
- Wrapped separately with PE shrink film



Suggested use

- MEDICAL
- BEAUTY
- HEALTH
- TATTOO
- HO.RE.CA
- CLEANING
- HOME USE



	REF	WEIGHT	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/CARTON	CARTONS/PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
Airlaid Examination Rolls	108.113.AIR	50 DS	50 cm	25 m	● PRINT	1	38 cm	66	12	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.115.AIR	50 DS	58 cm	25 m	● PRINT	1	38 cm	66	12		40 x 30 x 59 cm	80 x 120 x 255 cm
Airlaid & Pe Examination Rolls	108.123.AIR	50 DS	50 cm	25 m	● PRINT	2	38 cm	66	12	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.125.AIR	50 DS	58 cm	25 m	● PRINT	2	38 cm	66	12		40 x 30 x 59 cm	80 x 120 x 255 cm
Airlaid Industrial Rolls	108.151.AIR	500 gr	25 cm	40 m	● PRINT	1	38 cm	105	6	98	36 x 24 x 25 cm	80 x 120 x 190 cm
	108.153.AIR	1000 gr	29 cm	68 m	● PRINT	1	38 cm	179	6	63	48 x 32 x 29 cm	80 x 120 x 190 cm
	108.154.AIR	2000 gr	29 cm	138 m	● PRINT	1	38 cm	363	2	32	40 x 30 x 29 cm	80 x 120 x 255 cm

Custom width & length and perforation upon request.
For further information please contact us.

Ability to print other colors or designs
(e.g. logo) upon request.

All packaging dimensions may vary +/- 1% | All roll dimensions may vary +/- 2%
Palletization heights are suggested for transportation cost efficiency for your orders.

	REF	WEIGHT	HEIGHT	WIDTH	COLOURS	PLIES	PIECES	PACKS/CARTON	CARTONS/PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
Airlaid Wipes	108.101	50 DS	30 cm	33 cm	● PRINT	1	50	24	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.102	50 DS	30 cm	40 cm	● PRINT	1	50	24	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.103	50 DS	40 cm	40 cm	● PRINT	1	50	18	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.104	50 DS	40 cm	70 cm	● PRINT	1	50	12	24	40 x 30 x 69 cm	80 x 120 x 225 cm
	108.104.80	80 DS	40 cm	70 cm	● PRINT	1	50	6	24	40 x 30 x 69 cm	80 x 120 x 225 cm

SERIES	TYPES	COMPOSITION	ICONS	PLIES	LAMINATION	SOFTNESS	ABSORBENCY	RESILIENCE	PE	WATERPROOF / HYDROPHOBIC	COLOURS	SUGGESTED USE
ELITE STANDARD	Examination roll	100% Pure paper pulp & LDPE	  	2	Thermal	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM PLUS	Examination roll	100% Pure paper pulp & HDPE	  	2	Glue	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM STANDARD	Examination roll	100% Pure paper pulp & HDPE	  	2	Glue	●●●○○	●●●●○	●●●○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
ECO	Examination roll	100% Recycled Pulp & HDPE	  	2	Glue	●●○○○	●●●●○	●●○○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR STANDARD	Examination roll	100% Pure paper pulp	   	2	Glue	●●●○○	●●●●○	●●○○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR NORMAL	Examination roll	100% Pure paper pulp	   	1	Glue	●●●○○	●●●●○	●●●○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR BASIC	Examination roll	100% Recycled Pulp	   	2	Glue	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC STANDARD	Examination roll	100% Pure paper pulp	  	2	–	●●●○○	●●●●○	●●○○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC NORMAL	Examination roll	100% Pure paper pulp	  	1	–	●●●○○	●●●●○	●●●○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC BASIC	Examination roll	100% Recycled Pulp	   	2	–	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME EXTRA PLUS	Examination roll	100% Polypropylene	   	1	–	●●●●●	○○○○○	●●●●●	○	Hydrophobic	○●	MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME PLUS	Examination roll	100% Polypropylene	   	1	–	●●●●●	○○○○○	●●●●○	○	Hydrophobic	○●●●	MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME STANDARD	Examination roll	100% Polypropylene	   	1	–	●●●●●	○○○○○	●●●○○	○	Hydrophobic	○●● PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
PRESTIGE	Examination roll	100% Polypropylene & LDPE	   	2	Thermal	●●●●●	○○○○○	●●●●●	●	Hydrophylic / Waterproof (PE side)	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid & HDPE	  	2	Thermal	●●●●○	●●●●●	●●●●●	●	Hydrophylic / Waterproof (PE side)	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid	  	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Cleaning Wipes	100% Airlaid	  	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO, HO.RE.CA., CLEANING, HOME USE
GENTLE Alfa-Alfa	Cellulose wadding paper	100% Pure paper pulp	  	2	–	●●●●●	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Extra-Soft	Cellulose wadding paper	Recycled pulp	  	2	–	●●●●○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Alfa-Soft	Cellulose wadding paper	Recycled pulp	  	2	–	●●●○○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
DELICATE	Industrial paper roll	100% Pure paper pulp	  	2	–	●●●○○	●●●●○	●●●●○	○	Hydrophobic	○	HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE BASIC	Industrial paper roll	100% Recycled Pulp	  	2	–	●●●○○	●●●○○	●●●○○	○	Hydrophobic		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE	Industrial paper wipes	100% Pure paper pulp	 	2	–	●●●○○	●●●○○	●●●●○	○	Hydrophobic	○	HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING