

• **ECO**

Eco paper & PE examination rolls

- **2ply:** Recycled paper & HDPE
- Glue laminated
- Extra soft
- Highly resilient
- Waterproof
- Wet strength
- Perforated




Suggested use

- MEDICAL
- BEAUTY
- HEALTH
- TATTOO



eco-white

STANDARD	REF	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/ CARTON	CARTONS/ PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
	108.021.SE		50 m		2	38 cm	132	12			
	108.021.SEM80	39 cm	80 m	PRINT	2	38 cm	210	8	40	40 x 30 x 41 cm	80 x 120 x 220 cm
	108.021.SEM100		100 m		2	38 cm	264	6			
	108.023.SE		50 m		2	38 cm	132	12			
	108.023.SEM80	50 cm	80 m	PRINT	2	38 cm	210	8	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.023.SEM100		100 m		2	38 cm	264	6			
	108.025.SE		50 m		2	38 cm	132	12			
	108.025.SEM80	58 cm	80 m	PRINT	2	38 cm	210	8	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.025.SEM100		100 m		2	38 cm	264	6			
	108.027.SE		50 m		2	38 cm	132	12			
		68 cm	80	PRINT	2	38 cm	210	8	24	40 x 30 x 69 cm	80 x 120 x 225 cm
	108.027.SEM100		100 m		2	38 cm	264	6			



Custom width & length and perforation upon request.
For further information please contact us.

Ability to print other colors or designs
(e.g. logo) upon request.

All packaging dimensions may vary +/- 1% | All roll dimensions may vary +/- 2%
Palletization heights are suggested for transportation cost efficiency for your orders.

• **PREMIUM**

Paper & PE examination rolls

- **2ply:** 100% pure paper pulp & HDPE
- Glue laminated
- Super soft & highly resilient
- Waterproof
- Extra wet strength
- Perforated
- Wrapped separately with PE shrink film



Suggested use

- MEDICAL
- BEAUTY
- HEALTH
- TATTOO



STANDARD	REF	COLOURS	HEIGHT	LENGTH	PLIES	PERFORATED PER	PIECES	ROLLS/ CARTON	CARTONS/ PALLET	CARTONS DIMENSIONS (cm)	PALLET DIMENSIONS (cm)
	108.021	white	39 cm	50 m	2	38 cm	132	12	40	40 x 30 x 41	80 x 120 x 220
	108.023	green		50 m	2	38 cm	132	12	32		
	108.023.SAM80	light blue	50 cm	80 m	2	38 cm	210	8	32	40 x 30 x 51	80 x 120 x 219
	108.023.SAM100	pink		100 m	2	38 cm	264	6	32		
	108.025	purple		50 m	2	38 cm	132	12	32		
	108.025.SAM80	yellow	58 cm	80 m	2	38 cm	210	8	32	40 x 30 x 59	80 x 120 x 255
	108.025.SAM100	black		100 m	2	38 cm	264	6	32		
	108.027	kids		50 m	2	38 cm	132	12	24		
	108.027.SAM80	custom	68 cm	80 m	2	38 cm	210	8	24	40 x 30 x 69	80 x 120 x 225
	108.027.SAM100	PLUS		100 m	2	38 cm	264	6	24		



Custom width & length and perforation upon request.
For further information please contact us.

Ability to print other colors or designs
(e.g. logo) upon request.

All packaging dimensions may vary +/- 1% | All roll dimensions may vary +/- 2%
Palletization heights are suggested for transportation cost efficiency for your orders.

SERIES	TYPES	COMPOSITION	ICONS	PLIES	LAMINATION	SOFTNESS	ABSORBENCY	RESILIENCE	PE	WATERPROOF / HYDROPHOBIC	COLOURS	SUGGESTED USE
ELITE STANDARD	Examination roll	100% Pure paper pulp & LDPE	   	2	Thermal	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM PLUS	Examination roll	100% Pure paper pulp & HDPE	   	2	Glue	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM STANDARD	Examination roll	100% Pure paper pulp & HDPE	   	2	Glue	●●●○○	●●●●○	●●●○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
ECO	Examination roll	100% Recycled Pulp & HDPE	   	2	Glue	●●○○○	●●●●○	●●○○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR STANDARD	Examination roll	100% Pure paper pulp	    	2	Glue	●●●○○	●●●●○	●●○○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR NORMAL	Examination roll	100% Pure paper pulp	    	1	Glue	●●●○○	●●●●○	●●●○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR BASIC	Examination roll	100% Recycled Pulp	    	2	Glue	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC STANDARD	Examination roll	100% Pure paper pulp	   	2	–	●●●○○	●●●●○	●●○○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC NORMAL	Examination roll	100% Pure paper pulp	   	1	–	●●●○○	●●●●○	●●●○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC BASIC	Examination roll	100% Recycled Pulp	    	2	–	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME EXTRA PLUS	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●●●	○	Hydrophobic	○ ●	MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME PLUS	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●●○	○	Hydrophobic	○ ● ● ●	MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME STANDARD	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●○○	○	Hydrophobic	○ ● PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
PRESTIGE	Examination roll	100% Polypropylene & LDPE	    	2	Thermal	●●●●●	○○○○○	●●●●●	●	Hydrophylic / Waterproof (PE side)	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid & HDPE	   	2	Thermal	●●●●○	●●●●●	●●●●●	●	Hydrophylic / Waterproof (PE side)	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid	   	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Cleaning Wipes	100% Airlaid	   	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO, HO.RE.CA., CLEANING, HOME USE
GENTLE Alfa-Alfa	Cellulose wadding paper	100% Pure paper pulp	   	2	–	●●●●●	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Extra-Soft	Cellulose wadding paper	Recycled pulp	   	2	–	●●●●○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Alfa-Soft	Cellulose wadding paper	Recycled pulp	   	2	–	●●●○○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
DELICATE	Industrial paper roll	100% Pure paper pulp	  	2	–	●●●○○	●●●●○	●●●●○	○	Hydrophobic	○	HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE BASIC	Industrial paper roll	100% Recycled Pulp	  	2	–	●●●○○	●●●○○	●●●○○	○	Hydrophobic		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE	Industrial paper wipes	100% Pure paper pulp	 	2	–	●●●○○	●●●○○	●●●●○	○	Hydrophobic	○	HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING