

REGULAR

Embossed paper examination rolls

- **1 ply** or **2ply**: 100% pure paper pulp
- Glue laminated
- Embossed
- Soft
- Maximum absorbency
- Perforated
- Wrapped separately with PE shrink film



Suggested use

- MEDICAL
- BEAUTY
- HEALTH
- TATTOO



white



STANDARD NORMAL

	REF	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/ CARTON	CARTONS/ PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
STANDARD	108.011.SG		50 m		2	38 cm	132	12			
	108.011.SGM80	39 cm	80 m	PRINT	2	38 cm	210	8	40	40 x 30 x 41 cm	80 x 120 x 220 cm
	108.011.SGM100		100 m		2	38 cm	264	6			
	108.013.SG		50 m		2	38 cm	132	12			
	108.013.SGM.80	50 cm	80 m	PRINT	2	38 cm	210	8	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.013.SGM100		100 m		2	38 cm	264	6			
	108.015.SG		50 m		2	38 cm	132	12			
	108.015.SGM.80	58 cm	80 m	PRINT	2	38 cm	210	8	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.015.SGM100		100 m		2	38 cm	264	6			
	108.017.SG		50 m		2	38 cm	132	12			
	108.017.SGM.80	68 cm	80 m	PRINT	2	38 cm	210	8	24	40 x 30 x 69 cm	80 x 120 x 225 cm
	108.017.SGM100		100 m		2	38 cm	264	6			
NORMAL	108.011.SG1.NR		50 m		1	38 cm	132	12			
	108.011.SG1.NRM80	39 cm	80 m	PRINT	1	38 cm	210	8	40	40 x 30 x 41 cm	80 x 120 x 220 cm
	108.011.SG1.NRM100		100 m		1	38 cm	264	6			
	108.013.SG1.NR		50 m		1	38 cm	132	12			
	108.013.SG1.NRM80	50 cm	80 m	PRINT	1	38 cm	210	8	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.013.SG1.NRM100		100 m		1	38 cm	264	6			
	108.015.SG1.NR		50 m		1	38 cm	132	12			
	108.015.SG1.NRM80	58 cm	80 m	PRINT	1	38 cm	264	4	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.015.SG1.NRM100		100			38 cm	264	6			
	108.017.SG1.NR	68 cm	50 m	PRINT	1	38 cm	132	12	24	40 x 30 x 69 cm	80 x 120 x 225 cm

Custom width & length and perforation upon request.
For further information please contact us.

Ability to print other colors or designs
(e.g. logo) upon request.

All packaging dimensions may vary +/- 1% | All roll dimensions may vary +/- 2%
Palletization heights are suggested for transportation cost efficiency for your orders.

REGULAR

Embossed paper examination rolls

- **2ply**: 100% recycled paper pulp
- Glue laminated
- Embossed
- Soft
- Maximum absorbency
- Perforated
- Wrapped separately with PE shrink film



Suggested use

- MEDICAL
- BEAUTY
- HEALTH
- TATTOO



eco white



BASIC

	REF	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/ CARTON	CARTONS/ PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
BASIC	108.011.SGEM40		40 m					12			
	108.011.SGEM70	39 cm	70 m	eco white PRINT	2	38 cm	132	8	40	40 x 30 x 41 cm	80 x 120 x 220 cm
	108.011.SGEM90		90 m					6			
	108.013.SGEM40		40 m					12			
	108.013.SGEM70	50 cm	70 m	eco white PRINT	2	38 cm	132	8	32	40 x 30 x 51 cm	80 x 120 x 223 cm
	108.013.SGEM90		90 m					6			
	108.015.SGEM40		40 m					12			
	108.015.SGEM70	58 cm	70 m	eco white PRINT	2	38 cm	132	8	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.015.SGEM90		90 m					6			
	108.017.SGEM40		40 m					12			
	108.017.SGEM70	68 cm	70 m	eco white PRINT	2	38 cm	132	8	24	40 x 30 x 69 cm	80 x 120 x 225 cm
	108.017.SGEM90		90 m					6			

Custom width & length and perforation upon request.
For further information please contact us.

Ability to print other colors or designs
(e.g. logo) upon request.

All packaging dimensions may vary +/- 1% | All roll dimensions may vary +/- 2%
Palletization heights are suggested for transportation cost efficiency for your orders.

SERIES	TYPES	COMPOSITION	ICONS	PLIES	LAMINATION	SOFTNESS	ABSORBENCY	RESILIENCE	PE	WATERPROOF / HYDROPHOBIC	COLOURS	SUGGESTED USE
ELITE STANDARD	Examination roll	100% Pure paper pulp & LDPE	   	2	Thermal	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM PLUS	Examination roll	100% Pure paper pulp & HDPE	   	2	Glue	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM STANDARD	Examination roll	100% Pure paper pulp & HDPE	   	2	Glue	●●●○○	●●●●○	●●●○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
ECO	Examination roll	100% Recycled Pulp & HDPE	   	2	Glue	●●○○○	●●●●○	●●○○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR STANDARD	Examination roll	100% Pure paper pulp	    	2	Glue	●●●○○	●●●●○	●●○○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR NORMAL	Examination roll	100% Pure paper pulp	    	1	Glue	●●●○○	●●●●○	●●●○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR BASIC	Examination roll	100% Recycled Pulp	    	2	Glue	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC STANDARD	Examination roll	100% Pure paper pulp	   	2	–	●●●○○	●●●●○	●●○○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC NORMAL	Examination roll	100% Pure paper pulp	   	1	–	●●●○○	●●●●○	●●●○○	○	Hydrophobic	○	MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC BASIC	Examination roll	100% Recycled Pulp	    	2	–	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME EXTRA PLUS	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●●●	○	Hydrophobic	○ ●	MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME PLUS	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●●○	○	Hydrophobic	○ ● ● ●	MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME STANDARD	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●○○	○	Hydrophobic	○ ● PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
PRESTIGE	Examination roll	100% Polypropylene & LDPE	    	2	Thermal	●●●●●	○○○○○	●●●●●	●	Hydrophylic / Waterproof (PE side)	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid & HDPE	   	2	Thermal	●●●●○	●●●●●	●●●●●	●	Hydrophylic / Waterproof (PE side)	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid	   	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Cleaning Wipes	100% Airlaid	   	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic	○ PRINT	MEDICAL, HEALTH, BEAUTY, TATTOO, HO.RE.CA., CLEANING, HOME USE
GENTLE Alfa-Alfa	Cellulose wadding paper	100% Pure paper pulp	   	2	–	●●●●●	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Extra-Soft	Cellulose wadding paper	Recycled pulp	   	2	–	●●●●○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Alfa-Soft	Cellulose wadding paper	Recycled pulp	   	2	–	●●●○○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
DELICATE	Industrial paper roll	100% Pure paper pulp	  	2	–	●●●○○	●●●●○	●●●●○	○	Hydrophobic	○	HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE BASIC	Industrial paper roll	100% Recycled Pulp	  	2	–	●●●○○	●●●○○	●●●○○	○	Hydrophobic		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE	Industrial paper wipes	100% Pure paper pulp	 	2	–	●●●○○	●●●○○	●●●●○	○	Hydrophobic	○	HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING