

◦ **SUPREME**

Non-woven examination rolls

- **1 ply:** Non-woven (TNT) spunbond
- Extra soft
- Super resilient
- Hydrophobic
- Perforated
- Wrapped separately with PE shrink film



Suggested use



MEDICAL



BEAUTY



HEALTH



TATTOO



◦ **PRESTIGE**

Non-woven & PE examination rolls

- **2ply:** Non-woven (TNT) spunbond & LDPE
- Thermal laminated
- Extra soft
- Super resilient
- Waterproof
- Super wet strength
- Wrapped separately with PE shrink film



Suggested use



MEDICAL



BEAUTY



HEALTH



TATTOO



	REF	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/ CARTON	CARTONS/ PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
STANDARD	108.033 108.033	50 cm	70 m		1	38 cm	184	12	32	40 x 30 x 50 cm	80 x 120 x 219 cm
	108.035 108.035	58 cm			1	38 cm	184	12	32	40 x 30 x 60 cm	80 x 120 x 255 cm
	108.037 108.035	68 cm			1	38 cm	184	12	24	40 x 30 x 70 cm	80 x 120 x 225 cm
PLUS	108.033.SD 108.033	50 cm	70 m		1	38 cm	184	9	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.035.SD 108.035	58 cm			1	38 cm	184	9	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.037.SD 108.035	68 cm			1	38 cm	184	9	24	40 x 30 x 69 cm	80 x 120 x 225 cm
EXTRA PLUS	108.033.SD40 108.033	50 cm	50 m		1	38 cm	132	6	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.035.SD40 108.035	58 cm			1	38 cm	132	6	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.037.SD40 108.035	68 cm			1	38 cm	132	6	24	40 x 30 x 69 cm	80 x 120 x 225 cm

Custom width & length and perforation upon request.
For further information please contact us.

Ability to print other colors or designs
(e.g. logo) upon request.

All packaging dimensions may vary +/- 1% | All roll dimensions may vary +/- 2%
Palletization heights are suggested for transportation cost efficiency for your orders.

	REF	HEIGHT	LENGTH	COLOURS	PLIES	PERFORATED PER	PIECES	ROLLS/ CARTON	CARTONS/ PALLET	CARTONS DIMENSIONS	PALLET DIMENSIONS
STANDARD	108.043.ST	50 cm	50 m		2	38 cm	132	12	32	40 x 30 x 51 cm	80 x 120 x 219 cm
	108.045.ST	58 cm			2	38 cm	132	12	32	40 x 30 x 59 cm	80 x 120 x 255 cm
	108.047.ST	68 cm			2	38 cm	132	12	24	40 x 30 x 69 cm	80 x 120 x 225 cm

Custom width & length and perforation upon request.
For further information please contact us.

All packaging dimensions may vary +/- 1%
All roll dimensions may vary +/- 2%



SERIES	TYPES	COMPOSITION	ICONS	PLIES	LAMINATION	SOFTNESS	ABSORBENCY	RESILIENCE	PE	WATERPROOF / HYDROPHOBIC	COLOURS	SUGGESTED USE
ELITE STANDARD	Examination roll	100% Pure paper pulp & LDPE	   	2	Thermal	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM PLUS	Examination roll	100% Pure paper pulp & HDPE	   	2	Glue	●●●○○	●●●●○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
PREMIUM STANDARD	Examination roll	100% Pure paper pulp & HDPE	   	2	Glue	●●●○○	●●●●○	●●●○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
ECO	Examination roll	100% Recycled Pulp & HDPE	   	2	Glue	●●○○○	●●●●○	●●○○○	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR STANDARD	Examination roll	100% Pure paper pulp	    	2	Glue	●●●○○	●●●●○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR NORMAL	Examination roll	100% Pure paper pulp	    	1	Glue	●●●○○	●●●●○	●●●○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
REGULAR BASIC	Examination roll	100% Recycled Pulp	    	2	Glue	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC STANDARD	Examination roll	100% Pure paper pulp	   	2	–	●●●○○	●●●●○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC NORMAL	Examination roll	100% Pure paper pulp	   	1	–	●●●○○	●●●●○	●●●○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
CLASSIC BASIC	Examination roll	100% Recycled Pulp	    	2	–	●●●○○	●●●○○	●●○○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME EXTRA PLUS	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●●●	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME PLUS	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●●○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
SUPREME STANDARD	Examination roll	100% Polypropylene	    	1	–	●●●●●	○○○○○	●●●○○	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
PRESTIGE	Examination roll	100% Polypropylene & LDPE	    	2	Thermal	●●●●●	○○○○○	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid & HDPE	   	2	Thermal	●●●●○	●●●●●	●●●●●	●	Hydrophylic / Waterproof (PE side)		MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Examination roll	100% Airlaid	   	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO
HYPER	Cleaning Wipes	100% Airlaid	   	1	–	●●●●○	●●●●●	●●●●●	○	Hydrophobic		MEDICAL, HEALTH, BEAUTY, TATTOO, HO.RE.CA., CLEANING, HOME USE
GENTLE Alfa-Alfa	Cellulose wadding paper	100% Pure paper pulp	   	2	–	●●●●●	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Extra-Soft	Cellulose wadding paper	Recycled pulp	   	2	–	●●●●○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
GENTLE Alfa-Soft	Cellulose wadding paper	Recycled pulp	   	2	–	●●●○○	●●●○○	○○○○○	○	Hydrophobic		MEDICAL, HEALTH, CLEANING
DELICATE	Industrial paper roll	100% Pure paper pulp	  	2	–	●●●○○	●●●●○	●●●●○	○	Hydrophobic		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE BASIC	Industrial paper roll	100% Recycled Pulp	  	2	–	●●●○○	●●●○○	●●●○○	○	Hydrophobic		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING
DELICATE	Industrial paper wipes	100% Pure paper pulp	 	2	–	●●●○○	●●●○○	●●●●○	○	Hydrophobic		HEALTH, HO.RE.CA., HOME USE, BEAUTY, TATTOO, CLEANING